

Pre-eclampsia: From prevention to diagnosis of high-risk pregnancies

Telle Ukonaho


Maternal-Fetal Health

Reproductive Health Business, Revvity

revvity



40 000 MOMS AND 500 000 BABIES LOSE THEIR LIVES EVERY YEAR



PRE-ECLAMPSIA **AFFECTS** BOTH MOTHER AND THE CHILD **ALSO LATER IN LIFE**

40 000 MOMS AND 500 000 BABIES LOSE THEIR LIVES EVERY YEAR




10 million
PREGNANCIES
IS AFFECTED



20%
of PRETERM
BIRTHS are due to
hypertensive
disorders during
pregnancy



over 2.5
million PRETERM
BIRTHS are
caused by pre-
eclampsia each
year




60%
of MATERNAL
DEATHS due to
preeclampsia are
preventable

AVOIDING THIS CONDITION WOULD BRING SUBSTANTIAL IMPROVEMENTS TO MATERNAL AND FETAL HEALTH AND TO COSTS

Pre-eclampsia introduction

- Pre-eclampsia is one of the leading causes of maternal and perinatal morbidity and mortality
- Symptoms start after 20th week of gestation
- Pre-eclampsia is a multisystem disorder of pregnancy, clinically defined by hypertension and protein in the urine. In the absence of proteinuria, the finding of maternal organ dysfunction can serve as an indicator of pre-eclampsia
- As the only cure for pre-eclampsia after onset is delivery, pre-eclampsia often leads to premature delivery of the fetus
- Both the mother and her baby are at increased risk of serious health complications and even death


AVOIDING THIS CONDITION WOULD BRING SUBSTANTIAL IMPROVEMENTS TO MATERNAL AND FETAL HEALTH



40% of patients with new-onset hypertension or new-onset proteinuria will develop pre-eclampsia.

Barton & Sibai. Obstet Gynecol. 2008

MAKING A BIG DIFFERENCE IN PREGNANCY CARE



Professor Ranjit Akolekar,
Medway Maritime Hospital, UK.

"The success of the pregnancy is dependent on three factors; the baby, the placenta and the mother.

Every one of them has to be in absolutely good health to be able to have a healthy outcome for both the mother and the baby.

If we want to improve the health of the mothers and babies, we need to devote time, effort, resources and energy into screening for placental health, because that is far more likely to improve pregnancy outcomes than anything else "

MAKING A BIG DIFFERENCE IN PREGNANCY CARE

Placental insufficiency is present in virtually every major obstetrical syndrome"

ROMERO, AM J OBSTET GYNECOL. 2011.

BABY



PLACENTA

MOTHER

PLACENTAL INSUFFICIENCIES, especially preeclampsia leading to ADVERSE OUTCOMES of pregnancy

ARE MUCH MORE COMMON THAN ANEUPLOIDIES COMBINED.


BOTH MOTHER AND BABY ARE AFFECTED both short and long-term



- **Trisomy 21 = 1 in 600 (0.17%)**
- **Trisomy 18 = 1:1200 (0.08%)**
- **Trisomy 13 = 1:4200 (0.02%)**
- **Stillbirth = 1:250 (0.40%)**
- **Early PE = 1:200 (0.5%)**
- **Preterm PE = 1:66 (1.5%)**
- **Total PE = 1:40 (2.5%)**
- **Severe FGR = 1:33 (3.0%)**


MATERNAL AND NEONATAL MORTALITY

Global causes of maternal death



Pre-eclampsia, a hypertensive disorder, is one of the leading causes of maternal mortality



Causes of Neonatal Mortality



Premature birth is number one cause of neonatal mortality



Pre-eclampsia often leads to premature delivery of the fetus

National, regional, and global estimates of preterm birth in 2020, with trends from 2010 (Lancet, 2023)



Prematurity numbers are not that different between countries.
...and have not been decreasing.


Effects of being survivor of premature: Morbidity and costs



Frequency of death and major, intermediate and minor morbidity. Data are n (%).


Outcome	Delivery gestational age (weeks)															P for trend
	All (n=8334)	23 (n=43)	24 (n=114)	25 (n=124)	26 (n=169)	27 (n=159)	28 (n=196)	29 (n=213)	30 (n=262)	31 (n=312)	32 (n=451)	33 (n=639)	34 (n=1058)	35 (n=1477)	36 (n=3117)	
Death	119 (1.4)	19 (44.2)	36 (31.6)	15 (12.1)	19 (11.2)	13 (8.2)	4 (2.0)	4 (1.9)	4 (1.5)	3 (1.0)	1 (0.2)	1 (0.2)	0 (0.0)	0 (0.0)	0 (0.0)	<.001
Major morbidity*	657 (7.9)	19 (44.2)	60 (52.6)	68 (54.8)	88 (52.1)	64 (40.3)	43 (21.9)	48 (22.5)	36 (13.7)	22 (7.1)	39 (8.7)	27 (4.2)	46 (4.4)	42 (2.8)	55 (1.8)	<.001
Minor morbidity†	3136 (37.6)	4 (9.3)	18 (15.8)	39 (31.5)	59 (34.9)	77 (48.4)	144 (73.5)	147 (69.0)	206 (78.6)	255 (81.7)	344 (76.3)	406 (63.5)	540 (51.0)	402 (27.2)	495 (15.9)	<.001
Survival without any of the above morbidities	4422 (53.1)	1 (2.3)	0 (0.0)	2 (1.6)	3 (1.8)	5 (3.1)	5 (2.6)	14 (6.6)	16 (6.1)	32 (10.3)	67 (14.9)	205 (32.1)	472 (44.6)	1033 (69.9)	2567 (82.4)	

Infographic for parents on the mortality results among premature infants



Without efficient prediction & prevention,
prematurity is still the primary management
and outcome of preterm preeclampsia

Traditional Pre-eclampsia screening for all – PE management for symptomatic



Pre-eclampsia screening
NO SYMPTOMS



For all

Identify women with
high risk for
pre-eclampsia

**Screening based on
maternal/obstetric risk factors**

Recent studies: only ~40% are detected

SYMPTOMS of pre-eclampsia may appear
from wk 20 onwards



For symptomatic

Test (symptomatic) women with
suspected pre-eclampsia

Diagnosing the preeclamptic with standard clinical tests

Are not specific to Pre-eclampsia,

Don't provide a prognosis of PE

Don't provide a prediction in the time,

Do not inform about fetal well being

60%

Of maternal
deaths due to
pre-eclampsia
are
preventable



Majority of long-term effects of PREECLAMPSIA SURVIVORS

3-4x times the risk of
developing high blood
pressure

2x times to develop
heart disease

higher risk to develop
diabetes


...could be prevented




Majority of long-term effects of PREECLAMPSIA SURVIVORS

- Babies who survive can have
- **breathing issues, intestinal (digestive) problems, and bleeding in their brains.**
- Long-term problems may include **developmental delay** (not meeting the developmental milestones for his or her age), **CP, ADHD, Neurological disorders and lower performance in school.**

...could be prevented




Improved Pre-eclampsia screening for all – PE management for symptomatic



Pre-eclampsia screening
NO SYMPTOMS




For all
Identify women with high risk for pre-eclampsia




SYMPTOMS of pre-eclampsia may appear
from wk 20 onwards



For symptomatic
Test (symptomatic) women with suspected pre-eclampsia



PREECLAMPSIA CARE IMPACT



NUMBER 1 GOAL

Complications of preeclampsia (hypertensive disorder of pregnancy)

SHORT-TERM



Hypertension, coma, convulsions, kidney failure, liver failure, DEATH



Preterm birth, STILLBIRTH

LONG-TERM



Cardiovascular disease, kidney damage, Liver disease, stroke, diabetes mellitus




Delayed cognitive development, autism, death, cardiovascular disease, Stroke, Diabetes mellitus

Preeclampsia survivors – health affected for lifetime

Global causes of death

Our World in Data



Data source: IHME, Global Burden of Disease (2024)
OurWorldInData.org/causes-of-death | CC BY

How to improve pregnancy care for preventing pre-eclampsia?

1ST PRE-ECLAMPSIA SCREENING



Traditional screening method is simple but poor


ACOG/NICE guidelines pre-eclampsia screening is based on maternal characteristics and medical history

Moderate-risk factors
Nulliparity
Age >40 years
Inter-pregnancy interval >10 years
BMI at first visit >35 kg/m²
Family history of PE

High-risk factors
History of hypertensive disease in previous pregnancy
Chronic kidney disease
Autoimmune disease
Diabetes mellitus
Chronic hypertension
Multi fetal pregnancy



Simple and easy method to perform, but the
DETECTION RATE IS POOR





Screening based on maternal/obstetric risk factors

Recent studies: **only ~40% are detected**

The first step to better detection combined screening with PLGF

ACOG/NICE guidelines pre-eclampsia screening is based on maternal characteristics and medical history




COMBINED SCREENING






EFFICACY OF SCREENING AND ASPIRIN in reducing preterm pre-eclampsia



1. Rolnik DL et al. ASPREE trial; Aspirin versus Placebo in Pregnancies at High Risk for Preterm Preeclampsia. *N Engl J Med.* 2017;377(7):613-622. 2. Wright D et al. Aspirin for Evidence-Based Preeclampsia Prevention trial: influence of compliance on beneficial effect of aspirin in prevention of preterm preeclampsia. *Am J Obstet Gynecol.* 2017; 217:685.e1-5. 3. Poon LC et al. Aspirin for Evidence-Based Preeclampsia Prevention trial: effect of aspirin in prevention of preterm preeclampsia in subgroups of women according to their characteristics and medical and obstetrical history. *Am J Obstet Gynecol* 2017;217:585.e1-5. 4. Wright et al. Aspirin for Evidence-Based Preeclampsia Prevention trial: effect of aspirin on length of stay in the neonatal intensive care unit. *Am J Obstet Gynecol.* 2018; 218:612.e1-6

EFFECT OF PRE-ECLAMPSIA

on healthcare costs and infant survival



The earlier a baby is born, the higher the cost of care.

Table 1. Short-term costs of pre-eclampsia to the United States health care system.^[11]

Estimated unit and total health care cost for pre-eclampsia patients in the United States, by gestational age at birth (2012) using California Office of Statewide Health Planning and Development and commercial claims data

Costs	<28 wks (3604)	28-33 wks (23,624)	34-36 wks (41,856)	37 wks or longer (87, 596)	All (156, 680)
Maternal cost per birth	\$29,131	\$24,063	\$19,692	\$17,021	\$19,075
Infant cost per birth	\$282,570	\$59,803	\$11,112	\$6013	\$21,847
Combined cost per birth	\$311,701	\$83,866	\$30,804	\$23,035	\$40,922
Total health care cost	\$1.2 billion	\$2.0 billion	\$1.3 billion	\$2.0 billion	\$6.4 billion
Total cost because of infant cost, %	91%	71%	36%	26%	

The earlier a baby is born, the higher the risk of death or serious disability.

- In 2019 preterm birth and low birth weight accounted for about **17% of infant deaths** in US.
- Babies who survive can have **breathing issues, intestinal (digestive) problems, and bleeding in their brains.**
- Long-term problems may include **developmental delay** (not meeting the developmental milestones for his or her age), **CP, ADHD, Neurological disorders and lower performance in school.**

Birth	Death <5 y	Cerebral palsy	Impaired work capacity
23-27w	80%	9.1%	10.6%
28-30w	40%	6.0%	8.2%
31-33w	11%	1.9%	4.2%
34-36w	2.3%	0.3%	2.4%
≥37w	0.6%	0.1%	1.7%

867,692 live births Norway 1967-1983
[Moster et al. NEJM 2008;359:262](#)



ELSEVIER

Contents lists available at ScienceDirect

Pregnancy Hypertension: An International Journal of Women's Cardiovascular Health

journal homepage: www.elsevier.com/locate/preghy



Review article

The diagnostic accuracy of the NICE risk-stratification algorithm in predicting pre-eclampsia: A systematic review with *meta-analysis*

James S Morris^{a,*}, Siraj Abualnaja^a, Miriam Baumgarten^b

^a University of Cambridge School of Clinical Medicine, Addenbrooke's Hospital NHS Foundation Trust, Hills Road, Cambridge CB2 0QQ, United Kingdom

^b Department of Obstetrics and Gynaecology, Addenbrooke's Hospital NHS Foundation Trust, Hills Road, Cambridge CB2 0QQ, United Kingdom

- This systematic review of the NICE risk-stratification algorithm. Incorporating data from 17 eligible studies and 258,278 pregnancies.
- The NICE risk-stratification algorithm for preeclampsia performs remarkably poorly as a tool by which to allocate Aspirin prophylaxis.
 - 56.7% of preeclampsia cases missing
 - Costing £265000/year unnecessary prescriptions of Aspirin
- Should be replaced with an alternative tools like FMF tool.

Comparison of diagnostic accuracy of early screening for pre-eclampsia by NICE guidelines and a method combining maternal factors and biomarkers: results of SPREE

M. Y. TAN^{1,2}, D. WRIGHT³, A. SYNGELAKI¹, R. AKOLEKAR^{1,4}, S. CICERO⁵, D. JANGA⁶, M. SINGH⁷, E. GRECO⁸, A. WRIGHT³, K. MACLAGAN⁹, L. C. POON^{1,10#} and K. H. NICOLAIDES^{1,2#}

NIHR | National Institute
for Health Research

Journals Library



Efficacy and Mechanism Evaluation

Volume 7 • Issue 8 • November 2020


ISSN 2050-4365

Mini-combined test compared with NICE guidelines for early risk-assessment for pre-eclampsia: the SPREE diagnostic accuracy study

Liona C Poon, David Wright, Steve Thornton, Ranjit Akolekar,
Peter Brocklehurst and Kypros H Nicolaides



Showed low compliance to treatment if no proper screening and counselling is performed



- Compliance to aspirin (Do high risk women take their tablets?)
 - 23% with NHS method (history only)
 - 85% with FMF method (Combined)

It is important that the women know their risk and it is based on specific assessment with measurement results.

ASPIRIN – potential candidate for prevention

- Timing?

Low-dose Aspirin from 12 wks


Meta analysis on prophylactic aspirin

31 randomized studies, 32217 patients

• Preeclampsia 0.90 (95% CI 0.84-0.97)


Askie *et al*, Lancet 2007

≤ 16 wks
(n=764)



0.47 (0.34-0.65)
21.3% vs. 9.3%

> 16 wks
(n=10,584)



0.81 (0.63-1.03)
8.1% vs. 7.3%

0 .2 .4 .6 .8 1.0 1.2 1.4 1.6 1.8 2.0

27 studies (11,348 women)

Bujold *et al*, 2010

Aspirin ≤16 wks: early vs late PE

Preterm PE 0.11 (0.04-0.33)

Term PE 0.98 (0.42-2.33)

0.01 0.1 1 1.0

5 studies (n=556)

Roberge *et al*, 2012

Perinatal death

≤ 16 wks
(n=923)



0.28 (0.11-0.68)
0.85% vs. 4.2%

> 16 wks
(n=9,458)



0.90 (0.71-1.14)
2.6% vs. 2.9%

0 .2 .4 .6 .8 1.0 1.2 1.4

Roberge *et al*, 2012

Significant difference in impact on pre-term PE and perinatal death, dependent on when Aspirin prophylaxis was started.


ASPIRIN – potential candidate for prevention

- The dose?

A randomised crossover design study comparing the pharmacokinetics and pharmacodynamics of two single Doses of ORal Aspirin (75 mg v 150mg) in pregnant women at risk of pre-eclampsia (DORA): implications on assessing aspirin response and patient adherence to therapy.

Raya Vinogradov^{*†§c}, Oisín N. Kavanagh^{*#o}, Jeremy Palmer[#], Paul Murphy[§], Emma Curtis[§], Farhad Kamali^{§#} & Stephen Robson^{†§}.


The doses and the effect to thromboxane B2 levels from blood samples, in different timepoints post ingestion of aspirin



- 150 mg dose SA inhibits thromboxane more efficiently than the 75 mg dose.
- With a dose of 150 mg SA, plasma SA concentrations are consistently above baseline levels up to 12 hours after oral dosing
- Study was also able to demonstrate that thromboxane levels were fully recovered after a wash-out period of one week.

A Randomized Controlled Study Comparing the Efficacy of 75mg Versus 150mg Aspirin for the Prevention of Preeclampsia in High-Risk Pregnant Women

Nishi Sinha¹, Shruti Singh¹, Mukta Agarwal², Pramod K. Manjhi¹, Rajesh Kumar¹, Sunil Kumar Singh¹, Aakanksha Priya¹



- With women who are at high risk of developing preeclampsia, Aspirin 150 mg once a day at bedtime is more effective than Aspirin 75 mg in preventing preeclampsia with similar adverse outcomes (NICU admission, IUGR, neonatal death, still birth, eclampsia, HELLP syndrome, placental abruption and pulmonary edema).


Aspirin delays preeclampsia in women with chronic hypertension

Emmanuel Bujold ¹, Daniel L Rolnik ², Liona Poon ³, Argyro Syngelaki ⁴, David Wright ⁵, Kypros H Nicolaides ⁴

Affiliations + expand

PMID: 40287082 DOI: 10.1016/j.ajog.2025.04.046

FIGURE
Cumulative risk of PE among participants taking aspirin or placebo



Secondary analysis of the ASPRE data to better estimate the role of aspirin in women with chronic hypertension.

Aspirin use is associated with a shift in the cumulative risk of preterm PE to later gestational age (P=.051).

Findings suggest that, in women with chronic hypertension, aspirin shifts the distribution of delivery with PE to a later gestational age, in a similar fashion to that in high-risk women without chronic hypertension.

Women with chronic hypertension often have PE with delivery very early in gestation. Aspirin delays PE by several weeks and certainly reduces perinatal morbidity and neonatal hospitalizations.

Aspirin does not specifically reduce the number of preterm PE.

Consequently, women with chronic hypertension benefit from aspirin during pregnancy by reducing the earliest cases of PE.

Aspirin in prevention of preeclampsia: Universal aspirin – Side effects

Aspirin use during pregnancy and the risk of bleeding complications: a Swedish population-based cohort study. Hastie, R. et al. (2021). *American journal of obstetrics and gynecology*, 224(1), 95.e1–95.e12.

Complication	Aspirin (n=4,088)	No aspirin (n=309,536)	Odds ratio (95%CI)
Intrapartum hemorrhage	117 (2.9%)	4,695 (1.5%)	1.63 (1.30,2.05)
Postpartum hemorrhage	411 (10.2%)	24,036 (7.8%)	1.23 (1.08,1.39)
Neonatal intracranial hemorrhage	3 (0.07%)	17 (0.01%)	9.66 (1.88, 49.48)

Problems: Aspirin use defined as any tablet taken for any reason

Systematic review and meta-analysis of adverse events of low-dose aspirin and clopidogrel in randomized controlled trials.

McQuaid, K. R., & Laine, L. (2006). *The American journal of medicine*, 119(8), 624–638

Complication (n=75,000 non-pregnant)	Aspirin	No aspirin	RR (95%CI)
Major bleeding	26,673	26,712	1.71 (1.41, 2.08)
Major gastrointestinal bleeding	28,686	28,719	2.07 (1.61, 2.66)
Intracranial bleeding	27,671	27,712	1.65 (1.12, 2.44)

Absolute annual increase: Major bleeding 0.13%, Gastrointestinal bleeding 0.12%, Intracranial bleeding 0.03%, non-pregnant)

Utero Exposure to Aspirin and Risk of Asthma in Childhood.

Chu, S. et al (2016). *Epidemiology (Cambridge, Mass.)*, 27(5), 726–731.

Complication	Aspirin (n=11,215)	No aspirin (n=8,677)	aOR (95%CI)
Childhood asthma	552 (5%)	352 (4%)	1.3 (1.1, 1.6)

The risk of childhood asthma was increased by 30%

In utero drug exposure and hearing impairment in 2-year-old children A case-control study using the EFEMERIS database.

Foch, C., et al (2018). *International journal of pediatric otorhinolaryngology*, 113, 192–197.

Complication	Abnormal hearing (n=1,245)	Normal hearing (n=28,046)	aOR (95%CI)
In-utero ASA exposure	28 (2.2%)	396 (1.4%)	1.61 (1.09, 2.37)

The risk of hearing disorders was increased by 60%

Implementation of First-Trimester Screening and Prevention of Preeclampsia: A Stepped Wedge Cluster-Randomized Trial in Asia. Nguyen-Hoang et al; FORECAST Collaborators *Circulation*. 2024 Oct 15;150(16):1223-1235.

Complication	High risk with aspirin	High risk without aspirin
Postpartum hemorrhage	2.05%	1.44%
Fetal structural defects	0.72%	0.64%
Fetal chromosomal defects	0.1%	0
At least one serious adverse event	3.15%	2.08%

The risk of PE was screened

Aspirin versus Placebo in Pregnancies at High Risk for Preterm Preeclampsia.

Rolnik, D. L., et al (2017)*The New England journal of medicine*, 377(7), 613–622.

ASPRE study

Complication	Aspirin (n=798)	Placebo (n=822)	P-value
Vaginal bleeding	29 (3.6%)	21 (2.6%)	0.27
Nasal bleeding	16 (2.0%)	27 (3.3%)	0.15
Abruption without PE	5 (0.62%)	9 (1.09%)	0.31
Neonatal IVH ≥ G2	2 (0.25%)	1 (0.12%)	0.55

The risk of PE was screened

ORIGINAL RESEARCH

Placental growth factor before 11 weeks for screening of preterm preeclampsia: The PreMoM study

Rocío López Mármol^{1,2}  | José Alejandro Ávila Cabreja³ | Teresa de Haro Romero⁴  |
Catalina de Paco Matallana⁵  | Olga Ocón Hernández^{2,6}  | Otilia González-Vanegas^{2,6}  |
Pilar Carretero Lucena^{2,6} | María Paz Carrillo¹ | Juan Luis Delgado⁵  | Valeria Rolle⁷ |
Uzay Gormus⁸ | Liza Orahá⁸ | María M. Gil^{9,10}  | Francisca S. Molina^{2,6}

¹Department of Obstetrics and Gynecology, Hospital Universitario Virgen de las Nieves, Granada, Spain

²Instituto de Investigación Biosanitaria ibs.GRANADA, Granada, Spain

³Fundación Pública Andaluza para la Investigación Biosanitaria Andalucía Oriental, Granada, Spain

⁴Clinical Laboratory Management Unit, Hospital Universitario Clínico San Cecilio, Granada, Spain

⁵Department of Obstetrics and Gynecology, Hospital Universitario Virgen de la Arrixaca, Murcia, Spain

⁶Department of Obstetrics and Gynecology, Hospital Universitario San Cecilio, Granada, Spain

⁷Biostatistics and Clinical Research Unit, Hospital Universitario de Torrejón, Madrid, Spain

⁸Department of Clinical Laboratory, Revvity Omics Sweden, Sollentuna, Sweden

⁹Department of Obstetrics and Gynecology, Hospital Universitario de Torrejón, Madrid, Spain

¹⁰School of Medicine, Universidad Francisco de Vitoria, Pozuelo de Alarcón, Madrid, Spain

8-10+6 weeks vs 11-13+6 weeks screening performance

- The best discrimination (AUROC: 0.863; 95% CI, 0.754–0.971) and calibration were achieved by the model using PlGF between 11 and 13+6 weeks.
- PlGF measured before 11 weeks did not improve preterm PE screening performance

Preterm preeclampsia screening and prevention: a comprehensive approach to implementation in a real-world setting

2025 Jan 15;25(1):32. doi: 10.1186/s12884-025-07154-6.

The study confirms the feasibility of integrating comprehensive preeclampsia screening into real-world clinical practices

- The integration of preeclampsia screening **had a minimal effect on the time required** for aneuploidy screening, with results obtained within a rapid turnaround time.
- While the FMF algorithm provides better performance when all parameters (UtAPI, PlGF, and PAPP-A) are included,
- FMF algorithm can deliver a personalized risk score even if one or both placental markers are missing. Therefore, the protocol is applicable to lower-resource settings where not all parameters may be available.




Table 2 Time used for each component of the screening test and turnaround time of the screening test


Time Interval	Median (p5, p95)
History Collection (min)	2.0 (1.0, 4.25)
Blood Pressure Collection (min)	3.0 (2.0, 4.0)
Nuchal Translucency Scan (min)	25 (20,40)
Uterine artery Doppler measurement (min)	2.0 (1.0,2.5)
Time between blood sample collection and blood sample receipt (days)	1(1,2)
Time between sample receipt and results report (days)	2(1,3)

Is it Time for Precision Screening in Preeclampsia?

JACC Adv. 2025 Mar;4(3):101585. doi: 10.1016/j.jacadv.2025.101585.


Improved first-trimester preeclampsia screening would

1. reduce costs for patients and hospitals by
 - decreasing unnecessary hospitalizations and
 - identifying early-onset preeclampsia patients that need close monitoring during pregnancy.
2. It would also alleviate patient anxiety and provide additional information to help clinicians to risk-stratify patients.



Utilizing the proposed models, providers can capture a larger portion of high-risk patients in order to allocate aspirin appropriately, leading to improved morbidity, maternal outcomes, and future cardiovascular health.


Effective screening reduces health care costs and improves quality of life



Pre-eclampsia screening



A row of four colored squares, each containing an icon and text. From left to right: a blue square with a clipboard icon and 'RECORD MEDICAL HISTORY'; a teal square with a blood drop icon and 'BLOOD TEST PIGF 1-2-3'; a yellow square with a heart and pulse icon and 'MEASURE BLOOD PRESSURE'; and a pink square with an ultrasound icon and 'MEASURE ULTRASOUND'. Plus signs are placed between the squares.



How to enhance pre-eclampsia management?

2T/3T PRE-ECLAMPSIA SHORT-TERM PREDICTION &
AID IN DIAGNOSIS



Pre-eclampsia related sign or symptoms are not specific for pre-eclampsia



25%
of women with gestational hypertension will progress to pre-eclampsia

Standard clinical care tests are not specific to pre-eclampsia

The standard clinical diagnostic assessments in case of suspected pre-eclampsia contain comprehensive set of tests


BUT, tests:

- Are not specific to Pre-eclampsia
- Don't provide a prognosis of PE
- Don't provide a prediction in the time
- Lead to overdiagnosis and unnecessary hospital admission
- May lead to incorrect patient management
- Do not inform about fetal well being
- May lead to delay in diagnosis



For symptomatic women PLGF based testing improves pre-eclampsia management

During 2T and 3T
PLGF and sFlt-1 are both
predictive and diagnostic
for pre-eclampsia



Biomarker levels also correlate
with severity of disease

PLGF BASED TESTING in addition to standard clinical assessment **ENHANCE CLINICAL CARE**



For symptomatic women PLGF based testing – two alternatives

PIGF alone
Concentration
based cut-off

**sFlt-1/PIGF
ratio**



Both options are
equally recommended
for clinical use


PIGF,
as single marker,
has a comparable
performance
to sFlt-1/PIGF ratio

Based on the studies
there is no need
for two analytes,
PIGF alone
is sufficient

PLGF alone,
with concentration
based cut-offs,
could provide a
simpler solution

PLGF alone,
use of one kit,
could offer more
affordable
alternative





Increasing
implementation and
interest

Preeclampsia SCREENING/DIAGNOSING growth 2025

– supporting signals



Supporting data

Increasing amount of data showing the efficacy of screening and added value of PIGF and sFlt-1:

- FORECAST
- PREVAL
- PRECISE
- PRAECIS
- (RANSPRE, IMPACT

GUIDELINE

Supporting guidelines

International guidelines FIGO, ISSHP, ISUOG,

Local guidelines (screening): Germany, Austria, Switzerland, Denmark, Estonia, Canada, Brazil...Updating: Spain, Sweden

Local guidelines (management): Lithuania, UK, Finland, Germany, Austria, Switzerland, Spain, Canada, Brazil



Increasing implementation

National implementation studies in:

- Denmark
- France
- UK

New countries:

- USA

GRANTS:

- UNITAID
- Gates





Global strategic priorities

Organization with supporting strategies to accelerate reduction of maternal, newborn and child health

Hypertensive disorders of pregnancy being acknowledge as one of the leading causes of maternal mortality and morbidity

Increasing interest for preeclampsia screening



	2017	2018	2019	2020	2021	2022	2023	2024	2025
Screening	financed by EU 2017 ASPRE project		FIGO Poland		ISSHP		ISUOG Brazil	Australia New Zealand	
Both		Czech	Romania Bulgaria	Italy Spain Thailand	Estonia	Japan Canada		Germany Austria Switzerland	
Management	UK (NICE) Brazil	Malaysia ISUOG ESC	South Africa		Finland FIGO ISSHP	Sri Lanka	Australia New Zealand Lithuania Scotland	USA	Denmark


Increasing interest for preeclampsia screening



Supporting data



Increasing implementation



Supporting guidelines



Global strategic priorities

GRANTS and open calls:

Call for proposals: Improving access to lifesaving tools for prevention, diagnosis, and management of pre-eclampsia and maternal anemia



Date posted
21 November 2024

Call status
CLOSED

Organizational strategies to accelerate reduction of maternal, newborn and child health

- Hypertensive disorders of pregnancy being acknowledge as one of the leading causes of maternal mortality and morbidity (16%)

Reducing the Burden of Preeclampsia



Share this content

Initiative
Grand Challenges


Date Open
Feb 11, 2025, 10:00

Date Closed
Mar 25, 2025, 11:30

Mar 8, 2025

Many pregnancy-related complications going undetected and untreated - WHO

Haemorrhage – severe heavy bleeding – and hypertensive disorders like preeclampsia are the leading causes of maternal deaths globally, according to





Original Research

The Implementation of Preeclampsia Screening and Prevention (IMPRESS) Study

J.M. Johnson MD, FRCSC ^a, Jennifer D. Walsh MD, FRCSC ^a, Nanette B. Okun MD, FRCSC ^f, Amy Metcalfe PhD ^a, Melanie L. Pastuck BN ^a, Connor M. Maxey BSc ^a, Nancy Soliman MD, FRCSC ^a, Houman Mahallati MD, FRCPC ^b, Verena H. Kuret MD, FRCSC ^a, Shannon J. Dwinell MD, FRCSC ^a, Rati Chada MD, FRCSC ^a, Candace P. O'Quinn MD, FRCSC ^a, Jaime Schacher MD, FRCSC ^a, David A. Somerset DM, FRCSC ^a, Kimiko Paterson MD, FRCPC ^b, Ian R. Karim MD, FRCSC ^a

Ultrasound Obstet Gynecol 2023; 61: 682–690

Published online 5 May 2023 in Wiley Online Library (wileyonlinelibrary.com). DOI: 10.1002/ultr.26403

This is an open access article under the terms of the Creative Commons Attribution License, which permits use, distribution and reproduction in any medium, provided the original work is properly cited.

Pre-eclampsia screening in Denmark: a validation study

I. RIISHEDE^{1,2}, L. RODE^{2,3}, L. SPERLING^{4,5}, M. O. P. SANDAGER^{7,8}, H. SKOV^{7,8}, S. R. WAGNER⁹, P. N. C. A. N. C. K.

RESEARCH

Preterm preeclampsia and prevention: a cohort study to implementation in

Stefania Ronzoni
Tianhua Huang^{2,4}

Ultrasound Obstet Gynecol 2023
Published online in Wiley Online

Validation of a prediction model for first-trimester prediction of preeclampsia using cohort from PREVAL study

M. M. GIL^{1,2}, D. CUENCA-GÓMEZ^{1,2}, V. ROLLE^{1,3}, M. PERTEGAL^{4,5,6}, C. DÍAZ⁷, R. REVELLO⁸, B. ADIEGO⁹, M. MENDOZA¹⁰, F. S. MOLINA^{11,12}, B. SANTACRUZ^{1,2}, Z. ANSBACHER-FELDMAN¹³, H. MEIRI¹⁴, R. MARTIN-ALONSO^{1,2}, Y. LOUZOUN¹³ and C. DE PACO MATA LLANA^{4,5,6}



Clinical implementation of pre-eclampsia screening in the first trimester of pregnancy

Adalgisa Cordisco
Giulia Vannucci^b

^a Division of Prenatal Diagnosis,
^b Fetal Medicine Unit, Department

DOI: 10.1111/1471-0528.16361
www.bjog.org

Implementati

ester ia: a clinical

Increasing amount of data showing the efficacy of screening, added value of PIGF and cost savings.

On-going:


- RANSPRE (Fr)
- IMPACT (Swe)
- STARSHIP (UK)
- PEPPI (Fin)
- FNIH (US)




First Trimester Screening and A Stepped Wedge Cluster-

PhD^{*}, Angela S.T. Tai, MSc, Duy-Anh
MD, PhD ^{id}, Arihiro Shiozaki, MD, PhD,
MD, PhD ^{id}, Yali Hu, MD, Bin Li, MD ^{id}, Aditya Kusuma, MD, PhD ^{id},
Pengbuan Yapan, MD, Arundhati Gosavi, MD, PhD, Mayumi Kaneko, MD ^{id}, Suchaya
Luewan, MD ^{id}, Tung-Yao Chang, MD ^{id}, Noppadol Chaayasit, MD ^{id}, Tongta
Nanthakomon, MD, Huishu Liu, MD ^{id}, Steven W. Shaw, MD, PhD, Wing Cheong Leung,
MD ^{id}, Zaleha Abdullah Mahdy, MD, PhD ^{id}, Angela Aguilar, MD, MSc, Hillary H.Y.
Leung, MBBS, BSc, Nikki M.W. Lee, MD ^{id}, So Ling Lau, MD ^{id}, Isabella Y.M. Wah, MD
^{id}, Xiaohong Lu, PhD, Daljit S. Sahota, PhD ^{id}, Marc K.C. Chong, PhD ^{id}, and Liona
C. Poon, MD ^{id} FORECAST Collaborators

PE testing maturity level - Readiness to adopt



Pre-eclampsia screening for all – PE management for symptomatic



Pre-eclampsia screening




For all pregnant in 1T
Identify pregnant women
who has high risk for
pre-eclampsia

Symptoms of pre-eclampsia may appear
from wk 20 onwards




Test (symptomatic) women with
suspected pre-eclampsia later in pregnancy
(3-10% of pregnant)
Short-term prediction & aid in diagnosis

Pre-eclampsia screening for all – PE management for symptomatic



Preeclampsia affects both mother and the child also later in life

PREVENT




Rolnik et al. N Engl J Med 2017;377:613

BETTER OUTCOMES

Birth	Death <5 y	Cerebral palsy	Impaired work capacity
23-27w	80%	9.1%	10.6%
28-30w	40%	6.0%	8.2%
31-33w	11%	1.9%	4.2%
34-36w	2.3%	0.3%	2.4%
≥37w	0.6%	0.1%	1.7%

867,692 live births Norway 1967-1983
Moster et al. NEJM 2008;359:262



REDUCED COSTS

Cost of NICU Care

Gestational Age (completed weeks)	N	Mean (US \$)
24	486	222,563
25	678	233,538
26	756	207,637
27	900	178,080
28	1,091	146,121
29	1,226	115,801
30	1,556	92,882
31	1,995	68,446
32	2,799	46,117
33	4,719	30,145
34	14,541	10,535
35	25,077	6,007
36	44,922	3,444
37	92,421	2,027

Prophylactic use of aspirin results in

68% reduction

in the length of stay in NICU

MANAGE BETTER




RAPID & ACCURATE TRIAGE

THIS CONDITION WOULD BRING SUBSTANTIAL IMPROVEMENTS TO MATERNAL AND FETAL HEALTH AND COSTS

Time and commitment are needed in


BUILDING THE EVIDENCE




revvity

"Products may not be licensed in accordance with the laws in all countries, such as the United States and Canada. Please check with your local representative for availability."

Pre-eclampsia




SCREENING & PREVENTION - ALL WOMEN



PRE-ECLAMPSIA MANAGEMENT



Customer-Driven Innovation and Continuous Improvement




- Comprehensive
Multiple screening strategies are supported
- Effective data management with LifeCycle
Simplifies the data collection and workflow processes
- Easy to set up, learn and use
Simple and flexible user interface with performance monitoring dashboard

LifeCycle

Leading risk calculation software for prenatal screening

- Trusted system since 2000
- Over 800 users
- Supports all prenatal screening models
- Enables aneuploidy screening and pre-eclampsia screening



revvity

Products may not be licensed in accordance with the laws in all countries, such as the United States and Canada. More information: contact your representative.




Simple, flexible and reliable LifeCycle™ Prenatal Screening risk calculation software

- ✓ Complete data management solution for prenatal screening laboratories
- ✓ Simple and flexible user interface with a performance monitoring dashboard
- ✓ Comprehensive – Multiple screening strategies are supported
- ✓ Tailored to population – Optimized risk calculation to your population
- ✓ Accessible – Remote access for external users and connectivity to 3rd party applications
- ✓ Trusted solution over 20 years



Products may not be licensed in accordance with the laws in all countries, such as the United States and Canada. More information: contact your representative.





Innovations for
future

DIFFERENT NEEDS – DIFFERENT SOLUTIONS

TARGET CUSTOMERS

FROM CENTRALIZED
TO DECENTRALIZED



Central laboratories
& large hospitals



District hospitals



Private Obstetrics




Remote Healthcare workers



TARGET COUNTRIES

FROM LOW TO HIGH INCOME COUNTRIES
REMOTE TO NEAR HOSPITAL




Training
Infrastructure

Throughput
Complexity
Costs
Time-to-Result




TARGET PATIENTS


FROM SCREENING TO MANAGEMENT



These tests are typically sent to a hospital laboratory & **it can take up to 24 hours to receive a result**



Data from Europe





RONIA™ Platform

- **Quantitative** Placental Growth Factor (PIGF) results
- **Portable** reader (can run on batter)
- **No maintenance/service**
- **No need for laboratory facilities**
- **Whole blood** = no pre-treatment (centrifugation /laboratory skills) necessary. Serum also an option.
- **10µl sample volume**
- **No sample shipping:** immediate results = fast treatment, effective care
- **Short assay time (30 min), <2 min in the reader.**
- **Lot specific calibration** enabled with RFID



RONIA™ Platform

Result simply from whole blood finger prick



1. Add buffer to reaction tube



2. Collect finger prick sample



3. Add sample to the reaction tube - mix




4. Incubate 15 min



5. Add reaction to the cassette
6. Incubate another 15 min and measure

PlGF assay Total time ~ 30 min

RONIA PIGF Advanced kit workflow



revvity

Thank you!